

Clarius Ultrasound Scanners

Clarius Ultrasound Scanners are wireless and work with a mobile app that is compatible with most iOS and Android smart devices. Clarius scanners are designed to be carried around for quick exams and to guide procedures such as nerve blocks and targeted injections.

Clarius Scanner Options



	Basic	Premium
M-Mode	●	●
2 Rechargeable Batteries	●	●
1 Charger with Global AC Adapter	●	●
Access to Clarius Cloud Platform	●	●
Virtual Phased Array (Clarius C3/C7)	●	●
Virtual Linear (Clarius C3)	●	●
Color Doppler	●	●
Power Doppler	●	●
Pulsed-wave Doppler	●	●
Heart Rate Monitor (Clarius C3)	●	●
DICOM	●	●
Needle Enhancement (Clarius L7)	●	●
Premium Cloud Options	●	●

● These options are available for a fee.

Clarius Cloud Options

	Basic	Pro	Enterprise
Storage	2 GB/Scanner Up to 10 GB	100 GB	1 TB
Users	Unlimited	Unlimited	Unlimited
Reviews	50	Unlimited	Unlimited
PDF Reports	Unlimited	Unlimited	Unlimited
Act One™	1	50	Unlimited
Monthly Fees	Included	Contact us	Contact us

* Prices are quoted in US dollars and are applicable globally.

** RF data storage is available only on Clarius Cloud Enterprise.